

# ASAH1 Antibody, HRP conjugated



PACO52355

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes the sphingolipid ceramide into sphingosine and free fatty acid, .

**Gene ID:**

ASAH1

**Uniprot**

Q13510

**Synonyms:**

acid, ceramidase (AC) (ACDase) (acid, CDase) (EC 3.5.1.23) (Acylsphingosine deacylase) (N-acylsphingosine amidohydrolase) (Putative 32 kDa heart protein) (PHP32) [Cleaved into: acid, ceramidase subunit alpha; acid, ceramidase subunit beta], ASAH1, ASAH

**Immunogen:**

Recombinant Human acid, ceramidase protein (22-395AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A