

# CUL4A Antibody, FITC conjugated



PACO52360

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CUL4A

**Uniprot**

Q13619

**Synonyms:**

Cullin-4A (CUL-4A), CUL4A

**Immunogen:**

Recombinant Human Cullin-4A protein (1-61AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A