

FOXI1 Antibody, HRP conjugated



PACO52383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional activator required for the development of normal hearing, sense of balance and kidney function. Required for the expression of SLC26A4/PDS, JAG1 and COCH in a subset of epithelial cells and the development of the endolymphatic system in the inner ear. Also required for the expression of SLC4A1/AE1, SLC4A9/AE4, ATP6V1B1 and the differentiation of intercalated cells in the epithelium of distal renal tubules.

Gene ID:

FOXI1

Uniprot

Q12951

Synonyms:

Forkhead box protein I1 (Forkhead-related protein FKHL10) (Forkhead-related transcription factor 6) (FREAC-6) (Hepatocyte nuclear factor 3 forkhead homolog 3) (HFH-3) (HNF-3/fork-head homolog 3), FOXI1, FKHL10 FREAC6

Immunogen:

Recombinant Human Forkhead box protein I1 protein (183-343AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A