

EARS2 Antibody, Biotin conjugated



PACO52433

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the attachment of glutamate to tRNA(Glu) in a two-step reaction: glutamate is first activated by ATP to form Glu-AMP and then transferred to the acceptor end of tRNA(Glu).

Gene ID:

EARS2

Uniprot

Q5JPH6

Synonyms:

Probable glutamate--tRNA ligase, mitochondrial (EC 6.1.1.17) (Glutamyl-tRNA synthetase) (GluRS), EARS2, KIAA1970

Immunogen:

Recombinant Human Probable glutamate--tRNA ligase, mitochondrial protein (42-166AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A