

SLC3A1 Antibody, HRP conjugated



PACO52443

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the high-affinity, sodium-independent transport of cystine and neutral and dibasic amino acid, (system B(0, +)-like activity). May function as an activator of SLC7A9 and be involved in the high-affinity reabsorption of cystine in the kidney tubule.

Gene ID:

SLC3A1

Uniprot

Q07837

Synonyms:

Neutral and basic amino acid, transport protein rBAT (NBAT) (D2h) (Solute carrier family 3 member 1) (b(0, +)-type amino acid, transport protein), SLC3A1, RBAT

Immunogen:

Recombinant Human Neutral and basic amino acid, transport protein rBAT protein (1-68AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A