

# SLC3A1 Antibody, Biotin conjugated



PACO52445

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the high-affinity, sodium-independent transport of cystine and neutral and dibasic amino acid, (system B(0, +)-like activity). May function as an activator of SLC7A9 and be involved in the high-affinity reabsorption of cystine in the kidney tubule.

**Gene ID:**

SLC3A1

**Uniprot**

Q07837

**Synonyms:**

Neutral and basic amino acid, transport protein rBAT (NBAT) (D2h) (Solute carrier family 3 member 1) (b(0, +)-type amino acid, transport protein), SLC3A1, RBAT

**Immunogen:**

Recombinant Human Neutral and basic amino acid, transport protein rBAT protein (1-68AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A