

CYP20A1 Antibody, FITC conjugated



PACO52456

Product Information

Size:	Protein Background:
50ug	membrane.
Reactivity:	Gene ID:
Human	CYP20A1
Source:	Uniprot
Rabbit	Q6UW02
Isotype:	Synonyms:
IgG	Cytochrome P450 20A1 (EC 1.14), CYP20A1
Applications:	Immunogen:
ELISA	Recombinant Human Cytochrome P450 20A1 protein (91-390AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A