

PACO52490

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

Gene ID:

actbb

Uniprot

Q7ZVF9

Synonyms:

Actin, cytoplasmic 2 (Beta-actin-2) [Cleaved into: Actin, cytoplasmic 2, N-terminally processed], actbb, bactin2

Immunogen:

Recombinant Zebrafish Actin, cytoplasmic 2 protein (2-375AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A