

# SPRED1 Antibody, Biotin conjugated



PACO52505

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Tyrosine kinase substrate that inhibits growth-factor-mediated activation of MAP kinase. Negatively regulates hematopoiesis of bone marrow.

**Gene ID:**

SPRED1

**Uniprot**

Q7Z699

**Synonyms:**

Sprouty-related, EVH1 domain-containing protein 1 (Spred-1) (hSpred1), SPRED1

**Immunogen:**

Recombinant Human Sprouty-related, EVH1 domain-containing protein 1 protein (123-311AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A