

# FAM19A3 Antibody



PACO52506

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAM19A3

**Uniprot**

Q7Z5A8

**Synonyms:**

Protein FAM19A3 (Chemokine-like protein Tafa-3), FAM19A3, Tafa3

**Immunogen:**

Recombinant Human Protein FAM19A3 protein (31-169AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A