

PHOSPHO1 Antibody



PACO52538

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Phosphatase that has a high activity toward phosphoethanolamine (PEA) and phosphocholine (PCho). Involved in the generation of inorganic phosphate for bone mineralization.

Gene ID:

PHOSPHO1

Uniprot

Q8TCT1

Synonyms:

Phosphoethanolamine/phosphocholine phosphatase (EC 3.1.3.75), PHOSPHO1


Immunogen:

Recombinant Human Phosphoethanolamine/phosphocholine phosphatase protein (1-267AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO52538 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO52538 at dilution of 1:100.