

TBC1D31 Antibody



PACO52554

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

centrosome.

Gene ID:

TBC1D31

Uniprot

Q96DN5

Synonyms:

TBC1 domain family member 31 (WD repeat-containing protein 67), TBC1D31, WDR67


Immunogen:

Recombinant Human TBC1 domain family member 31 protein (662-1001AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry analysis of human testis tissue using PACO52554 at dilution of 1:100.



Immunocytochemistry analysis of human glioma using PACO52554 at dilution of 1:100.