

# NECTIN2 Antibody, HRP conjugated



PACO52563

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probable cell adhesion protein.

**Gene ID:**

NECTIN2

**Uniprot**

Q92692

**Synonyms:**

Nectin-2 (Herpes virus entry mediator B) (Herpesvirus entry mediator B) (HveB) (Nectin cell adhesion molecule 2) (Poliovirus receptor-related protein 2) (CD antigen CD112), NECTIN2, HVEB PRR2 PVRL2

**Immunogen:**

Recombinant Human Nectin-2 protein (32-360AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A