

# USP6NL Antibody, HRP conjugated



PACO52639

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts as a GTPase-activating protein for RAB5A and RAB43. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR. Involved in retrograde transport from the endocytic pathway to the Golgi apparatus. Involved in the transport of Shiga toxin from early and recycling endosomes to the trans-Golgi network. Required for structural integrity of the Golgi complex.

**Gene ID:**

USP6NL

**Uniprot**

Q92738

**Synonyms:**

USP6 N-terminal-like protein (Related to the N-terminus of tre) (RN-tre), USP6NL, KIAA0019

**Immunogen:**

Recombinant Human USP6 N-terminal-like protein (387-693AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A