

SERPINI1 Antibody, Biotin conjugated



PACO52733

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serine protease inhibitor that inhibits plasminogen activators and plasmin but not thrombin. May be involved in the formation or reorganization of synaptic connections as well as for synaptic plasticity in the adult nervous system. May protect neurons from cell damage by tissue-type plasminogen activator.

Gene ID:

SERPINI1

Uniprot

Q99574

Synonyms:

Neuroserpin (Peptidase inhibitor 12) (PI-12) (Serpini1), SERPINI1, PI12

Immunogen:

Recombinant Human Neuroserpin protein (191-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A