

KIF17 Antibody, Biotin conjugated



PACO52749

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transports vesicles containing N-methyl-D-aspartate (NMDA) receptor 2B along microtubules.

Gene ID:

KIF17

Uniprot

Q9P2E2

Synonyms:

Kinesin-like protein KIF17 (KIF3-related motor protein), KIF17, KIAA1405 KIF3X

Immunogen:

Recombinant Human Kinesin-like protein KIF17 protein (500-777AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A