

POLL Antibody, HRP conjugated



PACO52767

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Repair polymerase. Involved in base excision repair (BER) responsible for repair of lesions that give rise to abasic (AP) sites in DNA. Has both DNA polymerase and terminal transferase activities. Has a 5'-deoxyribose-5-phosphate lyase (dRP lyase) activity.

Gene ID:

POLL

Uniprot

Q9UGP5

Synonyms:

DNA polymerase lambda (Pol Lambda) (EC 2.7.7.7) (EC 4.2.99) (DNA polymerase beta-2) (Pol beta2) (DNA polymerase kappa), POLL

Immunogen:

Recombinant Human DNA polymerase lambda protein (16-211AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A