

FANCF Antibody, HRP conjugated



PACO52787

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. May be implicated in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability.

Gene ID:

FANCF

Uniprot

Q9NPI8

Synonyms:

Fanconi anemia group F protein (Protein FAF), FANCF

Immunogen:

Recombinant Human Fanconi anemia group F protein (295-341AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A