

# FANCF Antibody, FITC conjugated



PACO52788

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. May be implicated in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability.

**Gene ID:**

FANCF

**Uniprot**

Q9NPI8

**Synonyms:**

Fanconi anemia group F protein (Protein FAF), FANCF

**Immunogen:**

Recombinant Human Fanconi anemia group F protein (295-341AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A