

EIF4ENIF1 Antibody, FITC conjugated



PACO52792

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Nucleoplasmic shuttling protein. Mediates the nuclear import of EIF4E by a piggy-back mechanism.

Gene ID:

EIF4ENIF1

Uniprot

Q9NRA8

Synonyms:

Eukaryotic translation initiation factor 4E transporter (4E-T) (eIF4E transporter) (Eukaryotic translation initiation factor 4E nuclear import factor 1), EIF4ENIF1

Immunogen:

Recombinant Human Eukaryotic translation initiation factor 4E transporter protein (340-630AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A