

Acss2 Antibody, FITC conjugated



PACO52844

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Activates acetate so that it can be used for lipid synthesis or for energy generation.

Gene ID:

Acss2

Uniprot

Q9QXG4

Synonyms:

Acetyl-coenzyme A synthetase, cytoplasmic (EC 6.2.1.1) (Acetate--CoA ligase) (Acetyl-CoA synthetase) (ACS) (AceCS) (Acyl-CoA synthetase short-chain family member 2) (Acyl-activating enzyme), Acss2, Acas2 Acecs1

Immunogen:

Recombinant Mouse Acetyl-coenzyme A synthetase, cytoplasmic protein (1-701AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A