

USP18 Antibody, FITC conjugated



PACO52888

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Can efficiently cleave only ISG15 fusions including native ISG15 conjugates linked via isopeptide bonds. Necessary to maintain a critical cellular balance of ISG15-conjugated proteins in both healthy and stressed organisms.

Gene ID:

USP18

Uniprot

Q9UMW8

Synonyms:

Ubl carboxyl-terminal hydrolase 18 (EC 3.4.19) (43 kDa ISG15-specific protease) (hUBP43) (ISG15-specific-processing protease) (Ubl thioesterase 18), USP18, ISG43

Immunogen:

Recombinant Human Ubl carboxyl-terminal hydrolase 18 protein (18-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A