

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Rab effector involved in exocytosis. May act as scaffold protein. Plays a role in dendrite formation by melanocytes.

Gene ID:

RIMS2

Uniprot

Q9UQ26

Synonyms:

Regulating synaptic membrane exocytosis protein 2 (Rab-3-interacting molecule 2) (RIM 2) (Rab-3-interacting protein 3), RIMS2, KIAA0751 RAB3IP3 RIM2


Immunogen:

Recombinant Human Regulating synaptic membrane exocytosis protein 2 protein (667-943AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of PC-3 cells using PACO52902 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).