

# mib1 Antibody, FITC conjugated



PACO52916

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

mib1

**Uniprot**

H9GYU5

**Synonyms:**

E3 ubiquitin-protein ligase mib1 (EC 2.3.2.27) (Protein mind bomb) (RING-type E3 ubiquitin transferase mib1), mib1, mib

**Immunogen:**

Recombinant Zebrafish E3 ubiquitin-protein ligase mib1 protein (1-274AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A