

# DUX4 Antibody



PACO52982

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Involved in transcriptional regulation. May regulate microRNA (miRNA) expression.

**Gene ID:**

DUX4

**Uniprot**

Q9UBX2

**Synonyms:**

Double homeobox protein 4 (Double homeobox protein 10), DUX4, DUX10

**Immunogen:**


Recombinant Human Double homeobox protein 4 protein (217-402AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO52982 at dilution of 1: 100.