

# ABT1 Antibody, HRP conjugated



PACO52991

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Could be a novel TATA-binding protein (TBP) which can function as a basal transcription activator. Can act as a regulator of basal transcription for class II genes.

**Gene ID:**

ABT1

**Uniprot**

Q9ULW3

**Synonyms:**

Activator of basal transcription 1 (hABT1) (Basal transcriptional activator), ABT1

**Immunogen:**

Recombinant Human Activator of basal transcription 1 protein (197-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A