

ABT1 Antibody, Biotin conjugated



PACO52993

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Could be a novel TATA-binding protein (TBP) which can function as a basal transcription activator. Can act as a regulator of basal transcription for class II genes.

Gene ID:

ABT1

Uniprot

Q9ULW3

Synonyms:

Activator of basal transcription 1 (hABT1) (Basal transcriptional activator), ABT1

Immunogen:

Recombinant Human Activator of basal transcription 1 protein (197-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A