

DLX6 Antibody, HRP conjugated



PACO53007

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

transcription factor activity, sequence-specific DNA binding, nervous system development, skeletal system development.

Gene ID:

DLX6

Uniprot

P56179

Synonyms:

Homeobox protein DLX-6, DLX6

Immunogen:

Recombinant Human Homeobox protein DLX-6 protein (110-171AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A