

VNN1 Antibody, HRP conjugated



PACO53063

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Amidohydrolase that hydrolyzes specifically one of the carboamide linkages in D-pantetheine thus recycling pantothenic acid, (vitamin B5) and releasing cysteamine.

Gene ID:

VNN1

Uniprot

O95497

Synonyms:

Pantetheinase (EC 3.5.1.92) (Pantetheine hydrolase) (Tiff66) (Vascular non-inflammatory molecule 1) (Vanin-1), VNN1

Immunogen:

Recombinant Human Pantetheinase protein (97-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A