

JUND Antibody, HRP conjugated



PACO53075

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor binding AP-1 sites.

Gene ID:

JUND

Uniprot

P17535

Synonyms:

Transcription factor jun-D, JUND

Immunogen:

Recombinant Human Transcription factor jun-D protein (1-136AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A