

# JUND Antibody, Biotin conjugated



PACO53077

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor binding AP-1 sites.

**Gene ID:**

JUND

**Uniprot**

P17535

**Synonyms:**

Transcription factor jun-D, JUND

**Immunogen:**

Recombinant Human Transcription factor jun-D protein (1-136AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A