

I79_001071 Antibody



PACO53078

Product Information

Size:

50ug

Reactivity:

Cricetulus griseus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

I79_001071

Uniprot

G3GTT2

Synonyms:

Immunogen:

Recombinant Cricetulus griseus C-C motif chemokine protein (25-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A