

PACO53082

---

## Product Information

**Size:**

50ug

**Reactivity:**

Enterobacteria phage M13

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May initiate with G7P the virion concomitant assembly-budding process, by interacting with the packaging signal of the viral genome. The assembly-budding takes place at the host inner membrane. In turn, G7P and G9P are present at the end of the filamentous virion that emerges first from the bacterial host.

**Gene ID:**

IX

**Uniprot**

P69538

**Synonyms:**

Tail virion protein G9P (Coat protein C, polypeptide II) (G9P), IX

**Immunogen:**

Recombinant Enterobacteria phage M13 Tail virion protein G9P protein (1-32aa).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A