

PACO53090

Product Information

Size:

50ug

Reactivity:

Aspergillus giganteus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This protein inhibits the growth of a variety of fungal species.

Gene ID:

afp

Uniprot

P17737

Synonyms:

Antifungal protein, afp

Immunogen:

Recombinant Aspergillus giganteus Antifungal protein (44-94AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A