

hupB Antibody



PACO53094

Product Information

Size:

50ug

Reactivity:

Salmonella typhi

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Histone-like DNA-binding protein which is capable of wrapping DNA to stabilize it, and thus to prevent its denaturation under extreme environmental conditions.

Gene ID:

hupB

Uniprot

P0A1R9

Synonyms:

DNA-binding protein HU-beta (HU-1) (NS1), hupB


Immunogen:

Recombinant Salmonella typhi DNA-binding protein HU-β protein (1-90AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: hupB antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 30 kDa. Observed band size: 30 kDa.