

PACO53162

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May modify the assembly dynamics of microtubules, such that microtubules are slightly longer, but more dynamic.

Gene ID:

EML4

Uniprot

Q96HE3

Synonyms:

Echinoderm microtubule-associated protein-like 4 (EMAP-4) (Restrictedly overexpressed proliferation-associated protein) (Ropp 120), EML4, C2orf2 EMAPL4


Immunogen:

Recombinant Human EML4 protein (1-62AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: A549 whole cell lysate, HeLa whole cell lysate. All lanes: EML4 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 109, 103 kDa. Observed band size: 109 kDa.