

WSB1 Antibody, Biotin conjugated



PACO53205

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes type II iodothyronine deiodinase/DIO2. Confers constitutive instability to HIPK2 through proteasomal degradation.

Gene ID:

WSB1

Uniprot

Q9Y6I7

Synonyms:

WD repeat and SOCS box-containing protein 1 (WSB-1) (SOCS box-containing WD protein SWiP-1), WSB1, SWIP1

Immunogen:

Recombinant Human WD repeat and SOCS box-containing protein 1 protein (201-421AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A