

CAMLG Antibody, HRP conjugated



PACO53231

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Likely involved in the mobilization of calcium as a result of the TCR/CD3 complex interaction. Binds to cyclophilin B.

Gene ID:

CAMLG

Uniprot

P49069

Synonyms:

Calcium signal-modulating cyclophilin ligand (CAML), CAMLG, CAML

Immunogen:

Recombinant Human Calcium signal-modulating cyclophilin ligand protein (1-189AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A