

# CTSZ Antibody



PACO53246

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Exhibits carboxy-monopeptidase as well as carboxy-dipeptidase activity.

**Gene ID:**

CTSZ

**Uniprot**

Q9UBR2

**Synonyms:**

Cathepsin Z (EC 3.4.18.1) (Cathepsin P) (Cathepsin X), CTSZ

**Immunogen:**


Recombinant Human Cathepsin Z protein (127-287AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO53246 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO53246 at dilution of 1:100.