

CTSZ Antibody, Biotin conjugated



PACO53249

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Exhibits carboxy-monopeptidase as well as carboxy-dipeptidase activity.

Gene ID:

CTSZ

Uniprot

Q9UBR2

Synonyms:

Cathepsin Z (EC 3.4.18.1) (Cathepsin P) (Cathepsin X), CTSZ

Immunogen:

Recombinant Human Cathepsin Z protein (127-287AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A