

PACO53274

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Gene ID:

GJA5

Uniprot

P36382

Synonyms:

Gap junction alpha-5 protein (Connexin-40) (Cx40), GJA5


Immunogen:

Recombinant Human Gap junction alpha-5 protein (227-358AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO53274 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO53274 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).