

# HSD17B7 Antibody, HRP conjugated



PACO53287

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Responsible for the reduction of the keto group on the C-3 of sterols.

**Gene ID:**

HSD17B7

**Uniprot**

P56937

**Synonyms:**

3-keto-steroid reductase (EC 1.1.1.270) (17-beta-hydroxysteroid dehydrogenase 7) (17-beta-HSD 7) (Estradiol 17-beta-dehydrogenase 7) (EC 1.1.1.62) (Short chain dehydrogenase/reductase family 37C member 1), HSD17B7, SDR37C1

**Immunogen:**

Recombinant Human 3-keto-steroid reductase protein (128-198AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A