

KCNK1 Antibody, Biotin conjugated



PACO53301

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KCNK1

Uniprot

O00180

Synonyms:

Potassium channel subfamily K member 1 (Inward rectifying potassium channel protein TWIK-1) (Potassium channel K2P1) (Potassium channel KCNO1), KCNK1, HOHO1
KCNO1 TWIK1

Immunogen:

Recombinant Human Potassium channel subfamily K member 1 protein (268-336AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A