

MSMB Antibody, HRP conjugated



PACO53323

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular space, nucleus.

Gene ID:

MSMB

Uniprot

P08118

Synonyms:

Beta-microseminoprotein (Immunoglobulin-binding factor) (IGBF) (PN44) (Prostate secreted seminal plasma protein) (Prostate secretory protein of 94 amino acid,) (PSP-94) (PSP94) (Seminal plasma beta-inhibin), MSMB, PRSP

Immunogen:

Recombinant Human Beta-microseminoprotein (27-107AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A