

# NCR1 Antibody, FITC conjugated



PACO53328

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

**Gene ID:**

NCR1

**Uniprot**

O76036

**Synonyms:**

Natural cytotoxicity triggering receptor 1 (Lymphocyte antigen 94 homolog) (NK cell-activating receptor) (Natural killer cell p46-related protein) (NK-p46) (NKp46) (hNKp46) (CD antigen CD335), NCR1, LY94

**Immunogen:**

Recombinant Human Natural cytotoxicity triggering receptor 1 protein (121-230AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A