

SERINC1 Antibody, Biotin conjugated



PACO53361

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Enhances the incorporation of serine into phosphatidylserine and sphingolipids.

Gene ID:

SERINC1

Uniprot

Q9NRX5

Synonyms:

Serine incorporator 1 (Tumor differentially expressed protein 1-like) (Tumor differentially expressed protein 2), SERINC1, KIAA1253 TDE1L TDE2

Immunogen:

Recombinant Human Serine incorporator 1 protein (331-387AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A