

# SPINT2 Antibody, HRP conjugated



PACO53363

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Inhibitor of HGF activator. Also inhibits plasmin, plasma and tissue kallikrein, and factor XIa.

**Gene ID:**

SPINT2

**Uniprot**

O43291

**Synonyms:**

Kunitz-type protease inhibitor 2 (Hepatocyte growth factor activator inhibitor type 2) (HAI-2) (Placental bikunin), SPINT2, HAI2 KOP

**Immunogen:**

Recombinant Human Kunitz-type protease inhibitor 2 protein (71-194AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A