

TOM1 Antibody, Biotin conjugated



PACO53385

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in intracellular trafficking. Probable association with membranes.

Gene ID:

TOM1

Uniprot

O60784

Synonyms:

Target of Myb protein 1, TOM1

Immunogen:

Recombinant Human Target of Myb protein 1 protein (306-487AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A