

CRAT Antibody, Biotin conjugated



PACO53417

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Carnitine acetylase is specific for short chain fatty acid, . Carnitine acetylase seems to affect the flux through the pyruvate dehydrogenase complex. It may be involved as well in the transport of acetyl-CoA into mitochondria.

Gene ID:

CRAT

Uniprot

P43155

Synonyms:

Carnitine O-acetyltransferase (Carnitine acetylase) (EC 2.3.1.7) (Carnitine acetyltransferase) (CAT) (CrAT), CRAT, CAT1

Immunogen:

Recombinant Human Carnitine O-acetyltransferase protein (181-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A