

PACO53466

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

extracellular exosome, integral component of plasma membrane, NAD+ ADP-ribosyltransferase activity, protein ADP-ribosylation.

Gene ID:

ART3

Uniprot

Q13508

Synonyms:

Ecto-ADP-ribosyltransferase 3 (EC 2.4.2.31) (ADP-ribosyltransferase C2 and C3 toxin-like 3) (ARTC3) (Mono(ADP-ribosyl)transferase 3) (NAD(P)(+)-arginine ADP-ribosyltransferase 3), ART3, TMART


Immunogen:

Recombinant Human Ecto-ADP-ribosyltransferase 3 protein (279-352AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO53466 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO53466 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).